Departamento de
Física

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
2 / B-VAT II INSTRUCTION MANUAL


Power On and Off

After all cords and cables have been connected, turn the Processor Module POWER Switch to OFF and plug the Processor Module Power Cord into a properly grounded AC outlet. Then turn the Monitor ON/OFF, Processor Module POWER and Hand Controller ON/OFF switches to ON. Thereafter, use only the Hand Controller ON/OFF Switch to turn power on and off.

Note: For accurate results, the B-VAT II must be calibrated for patient refracting distance according to instructions on page19.

Hand Controller Key and Switch Summary

Note: Before using the Hand Controller, adjust display brightness and viewing angle according to instructions on page18.

Primary Key Functions

- [ON]...[SHIFT] Display characters or figures corresponding to the acuity printed on the key. Also changes the characters or figures.
- [CHANGE] Changes displayed characters or figures. Cancels SHIFT mode.
- [R-G] Illuminates the red-green display.
- [L-R] Change the current chart to the next larger or smaller size and changes the displayed characters and figures. Pressed after [CHANGE], inserts and moves the pointer. Pressed after [CONTRAST], changes contrast.
- [SINGLE] Alternately selects single-character lines or multiple-character lines.
- [LINE] Alternately selects a single-line or a multiple-line display.
- [EQU] Alternately selects lines of equal or decreasing size.
- [TICK] Displays the astigmatic clock dial.
- [SLOAN] Selects Sloan letters, Snellen letters, Landolt rings or numbers.
- [CHILD] Selects children's symbols, tumbling E's, OC's or HOTV.
- [ACC.1] [ACC.2] [ACC.3] Alternately turns accessory outlets on or off.
- [SHIFT] Changes the key functions to those printed in yellow above each key. Cancelled by [CHANGE].
- [BLANK] Displays a blank screen on the Monitor.
- [FIXATE] Illuminates the Fixation Light.

Hand Controller.

*See Appendix C for information on other models.
Rev.1192.0

Quick Guide / 3

Alternate [SHIFT] Key Functions

- Printed in yellow above each key, these functions are activated by pressing [SHIFT]. (Press [CHANGE] to cancel the SHIFT mode.)
- [INVERSE] Alternately selects white characters on black or black characters on white.
 - [<=>] Displays the same characters one size larger.
 - [<=>] Displays the same characters one size smaller.
 - [SUM Rx] Shifts to the "Adding Prescription" mode.
 - [A] Not used; available for future options.
 - [B] Not used; available for future options.
 - [C] Not used; available for future options.
 - [TRANSP] Enters Menu mode. In Sum Rx, transposes a prescription from positive to negative cylinder value and vice versa.
 - [ASTIG] Displays the astigmatic "T".
 - [GRATING] Selects the grating optotype.
 - [CONTRST] Permits display contrast to be changed by pressing <=> or <=>.
 - [ALT DISP] Used with optional BVS Binocular Vision Testing System
 - [ANIMATE] Displays an animated fixation target.

Prescription-Addition Mode ([SHIFT] [SUM Rx]) Key Functions

- The following key functions are printed in red above each acuity key, and in white below the [R-G], [L-R], and [CHILD] keys. These functions are activated by pressing [SHIFT] and then [SUM Rx] once.
- The primary key functions are restored by repeatedly pressing [SHIFT] until the prescriptions are cleared from the display.
- [d]...[s] Enter prescription numeric values.
 - [+/-] Changes the sign of the value being entered.
 - [SPH] Enters the keyed value into the sphere position of the display.
 - [CYL] Enters the keyed value into the cylinder position of the display.
 - [AXIS] Enters the keyed value into the axis position of the display.
 - [C] Alternately displays the first and second prescriptions.

SPECIFICATIONS

Processor Module (Mentor Part No. 22-4810)

Dimensions: 14-1/4 L x 14-1/2 W x 4-1/4 H inches (36 x 37 x 11 cm)
 Weight: 13 lb. (6 kg.)
 Power requirements: 120 volts, 50/60 Hz, 100 watts (processor module alone)
 Accessory outlets: 350 watts max. total for all ACC outlets combined
 Maximum system power consumption: 600 watts including Monitor/Accessories

Monitor (Mentor Part No. 22-4805 or 22-4819)

Manufacturer's model no.: Electrohome Ltd. model V38-V30ILJ-AP
 Dimensions: 13-1/16 L x 12-1/2 W x 11-7/8 H inches (33.2 x 31.1 x 30.2 cm)
 CRT diagonal measurement: 11 in.
 Black Level Stability: 2% max. change from 10% to 90% APL 20 ft-L
 Picture Size Stability: 2% max. change 0 to 30 ft-L
 Geometric Distortion: 5% max.
 Linearity: 5% max.
 Weight: 25 lb. (11.3 kg.)
 Power requirements: 120/240 volts, 50/60Hz, 60 watts

Hand Controller (Mentor Part No. 22-4806)

Dimensions: 9 L x 2-3/4 W x 1-1/8 H inches (23 x 7 x 2.8 cm)
 Weight: 13 oz. (0.4 kg.)
 Controls: 27 illuminated function keys, one On/Off switch
 Display: Liquid crystal, illuminated, three lines, 10 characters/line
 LCD viewing angle: 15 - 90°

Rev. 1192.0

Appendix A

(OPTIONAL FEATURE - Available only on B-VAT II-SG Model)

CONTRAST SENSITIVITY TESTING WITH SINUSOIDAL GRATINGS

CONTRAST LEVELS AND GRATING SPATIAL FREQUENCIES

When equipped with the optional sinusoidal gratings, the B-VAT II can display any of the following 20 contrast levels and 16 grating spatial frequencies for use in contrast sensitivity testing:

Available Contrast Levels		Grating Spatial Frequencies	
Percent Contrast ⁽¹⁾	Equivalent Contrast Sensitivity ⁽²⁾	Cycles per Degree of Visual Angle	Equivalent Snellen Ratio
98	1.0	1.5	20/400
80	1.3	2.0	20/300
63	1.6	2.4	20/250
50	2.0	3.0	20/200
40	2.5	3.7	20/160
32	3.1	4.8	20/125
25	4.0	6.0	20/100
20	5.0	7.5	20/80
16	6.3	9.4	20/70
10	10	11.9	20/60
6.3	16	12.0	20/50
4.0	25	15.0	20/40
2.5	40	20.0	20/30
1.6	63	24.0	20/25
.63	160	30.0	20/20
.40	250	40.0	20/15
.25	400		
.16	630		
.10	1000		

⁽¹⁾ With the sinusoidal option installed and the gratings optotype displayed, the B-VAT II calculates percent contrast as follows:

$$\text{Percent Contrast} = 100 \times (\text{Max} - \text{Min}) / (\text{Max} + \text{Min})$$

These contrast levels differ from those listed under "Low Contrast Testing" on page 30 and they are calculated with a different formula. With the sinusoidal option installed, however, B-VAT II continues to calculate percent contrast according to the formula on page 30 whenever optotypes other than gratings are displayed.

⁽²⁾ Contrast Sensitivity = 100/(Percent Contrast)

Rev.0189.1

CHANGING CONTRAST AND SPATIAL FREQUENCY OF SINUSOIDAL GRATINGS

- Press **[SHIFT]**, then **[GRATING]** to enter the sinusoidal grating mode. When the pattern has reached the desired contrast, the Hand Controller will beep twice. When testing near threshold levels, this signal helps the examiner and the patient to know that the pattern has reached its desired contrast. The sound can be turned off (or back on) by pressing **[SHIFT]** **[C]**.

- Press **[SHIFT]**, then **[CONTRAST]** to enter the contrast adjustment mode. The monitor will blank to the average screen luminance during this adjustment.
- Press **[LEFT]** or **[RIGHT]** to change contrast one step at a time to the desired level. The percentage contrast will be shown on the Hand Controller display. Press **[CHANGE]** to bring back the display.

Note: When the display is at maximum contrast, pressing **[LEFT]** will change the contrast to the minimum value. At minimum contrast **[RIGHT]** changes the contrast to 98%.

In order to use the arrow keys for changing spatial frequency after contrast has been set to the desired level, press any acuity key or **[CHANGE]**. Spatial frequency can now be changed via the acuity keys or **[LEFT]** and **[RIGHT]** to determine the patient's best vision at the chosen contrast. The same contrast will be maintained until steps 2 and 3 are repeated to change contrast, or until full contrast is restored as described below.

RESTORING MAXIMUM CONTRAST

Return to a maximum contrast display by using any of these procedures:

- Repeat steps 2 and 3 in "Changing Contrast Levels" above to change the maximum contrast level.
- Press **[MAX]** twice.
- Turn the Hand Controller ON/OFF switch to OFF and back to ON.

Rev.0989.1

CONTROLLING THE TIMING PARAMETERS

Sinusoidal grating patterns are displayed on the patient's monitor with an edge gradient and with timing parameters in order to reduce certain characteristics that may bias the results of Contrast Sensitivity Testing.

When the timing parameters are in use, the monitor will present the gratings by increasing the contrast from zero percent contrast to the selected contrast, over a preset RISE TIME. The pattern will remain at the selected contrast level for a preset ON TIME and then gradually decrease to zero percent contrast over a preset FALL TIME. Before the next selected pattern is presented, a minimum OFF TIME must elapse before the new pattern begins its RISE TIME.

When the edge gradient is in use, the grating will increase from zero percent contrast to the selected contrast over a distance equal to about 0.2 degrees of visual angle. The total width of the illuminated pattern is 2.1 degrees.

The values for the timing parameters are as follows:

RISE TIME	0 to 9.75 seconds, adjustable in 0.25 second increments.
ON TIME	0.1 to 9 seconds, adjustable in 2 second increments down to 1 second and the following two steps are 0.20 and 0.1 second respectively. Invoking greater than 9 seconds will cause the gratings to remain on.
FALL TIME	0 to 9.75 seconds, adjustable in 0.25 second increments. "FALL TIME" does not apply when the grating "ON TIME" is indefinite.
OFF TIME	0 to 8 seconds, adjustable in 2 second increments. This is the minimum amount of time between grating patterns. "OFF TIME" does not apply when the grating "ON TIME" is indefinite.

The B-VAT II-SG is shipped from the factory with the timing parameters preset to the following values:

RISE TIME	0 seconds
ON TIME	Remains on
FALL TIME	0 seconds
OFF TIME	0 seconds

Rev.0989.1

The gradient timing parameters may be modified by using the **Parameter Set-up menu** (see page 17).

When modifying the timing parameters, the name of the parameter being adjusted, and its current value, are presented in the Hand Controller display. Press **CHANGE** to select the respective parameter then pressing the **UP** arrow or **DOWN** arrow will increment or decrement the current value of the parameter.

If desired, press **[.]** (decimal point) to store the current set of parameters as the default parameters to be in effect whenever the unit is powered up. Using **"SET DEFAULTS"** cancels all modified timing parameters and returns them to the preset factory values listed above.


While in gratings mode, **[SHIFT] [B]** toggles the edge gradient on and off.

NOTE: The screen luminance remains constant throughout the test at 85 candelas/m²

FURTHER INFORMATION

More information on contrast sensitivity testing is available from Mentor Customer Service. Call toll-free in the continental U.S. 800-992-7557.

Rev.0491.1

Rev.1192.0

Appendix E

(OPTIONAL FEATURE - Available only on B-VAT™ II-SG Model)

JAC™ JOYSTICK AUTOMATED CONTRAST SYSTEM (22-4871)

Introduction

The JAC system was developed to perform contrast sensitivity testing in a fast, automatic but flexible manner, with a printout of the patient's contrast sensitivity curve.

The system consists of two EPROMS, a Joystick with templates, an Encoder Box (to connect the Joystick and Hand Controller with the Processor Module) and a Printer for hard copy of the contrast sensitivity curve. This system may only be used with the B-VAT™ II-SG Video Acuity Tester.

Unique to the software in the system is the ability to program from 1 to 16 points to be tested on any of 31 contrast levels which range from 0.10% to 98% and change in equal 0.1 log unit steps. This represents an even greater contrast range than the 20 levels offered in the standard B-VAT II-SG model.

Unlike other systems, the JAC system enables testing with all nine of the optotypes found in the SG program. In addition to testing with sinusoidal gratings, the use of letters and other optotypes is possible with the new INCORRECT/CORRECT program. The program enhances accuracy while reducing the odds for correctly guessing a target. This testing format is also helpful where patients may be unable to use the Joystick. The technician determines whether a patient is correct based on his/her verbal response. Flexibility and versatility are built in to the system.

The JAC system also offers the ability to pre-program testing parameters. Timing parameters, the staircase method, and number of reversals can all be selected through the Hand Controller and a new software menu.

Using the JAC System

PERFORMING A FIVE POINT TEST

A five point test is currently programmed in the software. Changes can be made by using the menu. The five default starting points are as follows:

TEST	VA	CPD	CONTRAST
1	20/400*	1.5	6.3%
2	20/200	3.0	2.5%
3	20/100	6.0	2.5%
4	20/50	12.0	8.0%
5	20/30	20.0	25.0%

* 20/300 for non-gratings

TO PERFORM A 5 POINT TEST

Select optotype then:
Press

[SHIFT] [SUM Rx]

LCD reads:

★ Indicates Joystick must now be used for testing

▲ Indicates allowable of directions which can be used with the Joystick

Use appropriate template on the Joystick to facilitate patient response. Patient uses the Joystick to respond by moving the Joystick towards the graphic image which matches what appears on the screen.


Notes:

- Patient must respond by guessing even when unsure of target.
- Patient moves joystick to identify the target.
- System automatically selects next target based on response.
- System automatically cycles through all 5 points.
- Testing is over when a long "chirp" is heard.

Rev. 1192.0

SHIFT [C] 1st Press: Summary of thresholds

SHIFT [C] 2nd Press: Prints a graph of thresholds

A threshold is determined by averaging the log of the lowest contrast levels that are correct for the respective points. The results are rounded up to the next more visible step.

Rev. 0491.1

Select optotype desired.

Select desired acuity level to be tested (i.e. 20/70).

Select desired contrast at which to begin test.

Press **SHIFT** **ACC 2** to begin test.

Note: To speed up test enter a lower contrast level to be tested as the starting point. System will automatically begin the test from initial point selected.

If the Joystick has been selected to operate in the normal mode (not INCORRECT/CORRECT mode), when letters are presented, the JAC system will temporarily revert to the INCORRECT/CORRECT mode for the test.

CORRECT/INCORRECT template can now be used (Tubling E or HOTV can be used if desired) with Snellen Letter Test. Technician must now move the Joystick based on patient's response and information available on the LCD readout.

Technician moves Joystick to INCORRECT or CORRECT based on patient's verbal response.

The following single key strokes can be used.

The following single key strokes can be used.

[SUM Rx] Begin a Multi-Point test

[CLOCK] Begin a One-Point test

[C] 1st Press: Summary of thresholds
2nd Press: Prints a graph of thresholds

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

Apêndice - Manual do instrumento

- (1) EYE : Selection of an eye to be examined.
- (2) STAND.BY : It sets the screening level.
- (3) SCREEN : It selects the screening test.
- (4) ASSESS : It discriminates scotoma and make measurement of threshold value.
- (5) START/STOP : Press it to start examination. If it is pressed while examination is in process, the examination will be cancelled.
- (6) PAUSE : Press it to pause examination. If it is pressed not during examination, the central fovea lights with LED.
- (7) PRINT : Press it to print out the results while it is flickering.
- (8) REACTION LED : Green LED shows lighting of the targets.
Red LED shows response of a patient.
- (9) EXPOSURE : It shows the number of presentation (lighting) times of the targets to the blind spot.
- (10) RESPONSE : It shows the number of response times against the presented targets to the blind spot.
- (11) SCR.LEVEL : It shows the screening level.
- (12) THRESHOLD : It shows a threshold value and luminosity.
- (13) It increases or decreases the values of SCR. LEVEL and THRESHOLD.

Target presentation time select switch (0.2, 0.4 or 0.8 sec.)

Target presentation interval select switch (1, 2 or 3 sec.)

Main switch

Automatic print switch


When it is on, the final examination result can be printed out automatically even if an operator fails to press the PRINT button, and then the instrument gets "Stand-by" condition for the next patient.

Four (4) points fixation switch

When it is on, four fixation targets light up.


8. Preparation

- (1) Have a patient sit relaxedly in front of the instrument.
- (2) Cover an unexamined eye with an eyepatch, etc. (Examine a healthy or better eye first.)
- (3) Have the patient gaze at a fixation point with both eyes.
- (4) Have the patient hold a response push-button switch with his thumb on a push-button.
- (5) When a patient presses the response push-button correctly, a small electronic "click" sound is heard.
- (6) If a patient keeps on pressing the response push-button, the examination is interrupted and when he releases it, the examination starts again.
- (7) Set a patient's chin on a chin-rest and his forehead on a forehead rest properly. Adjust the eye height with a chin-rest vertical adjusting knob to be on a level with a white line. Make sure that a patient is in a relaxed posture.


When a correction lens is required, insert a lens holder and put the correction lens on it

- (8) Having prepared for the examination, present a target on a central fovea so that a patient can know the targets.


- (9) Tell a patient to concentrate on the fixation point and start examination.


- 22 -

- 25 -


- 27 -

- 26 -

11. Explanation on recording paper

The circle lit up with LED should be marked as a scotoma. A relative defect is shown with a flickering light. After measuring a threshold value, write down the value in the circle.

15. A sensitivity conversion table of ATS-95

The following is the conversion table of ATS-95 (dB, the left column) and Goldmann, Octopus and luminosity (asb). As the step of this instrument is 2dB, the actually used values are underlined. The column of the Goldmann shows the combination of the filters, 0 dB of Octopus is 1000 asb.

dB	asb	Goldmann	Octopus (dB)
0	40.0	III/4-a	4
1	31.5	3-e	5
2	25.0	d	6
3	20.0	b	7
4	16.0	a	8
5	12.5	2-e	9
6	10.0	d	10
7	8.0	c	11
8	6.3	b	12
9	5.0	a	13
10	4.0	1-e	14
11	3.15	d	15
12	2.5	c	16
13	2.0	b	17
14	1.6	a	18
15	1.25	0.01x4-e	19
16	1.0	d	20
17	0.8	c	21
18	0.63	b	22
19	0.5	a	23
20	0.4	3-e	24
21	0.315	d	25
22	0.25	c	26
23	0.2	b	27
24	0.16	a	28
25	0.125	2-e	29
26	0.1	d	30
27	0.08	c	31
28	0.063	b	32
29	0.05	a	33
30	0.04	1-e	34
31	0.0315	d	35
32	0.025	c	36
33	0.02	b	37
34	0.016	a	38

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Referências

http://www.psypress.co.uk/mather/resources/swf/Demo8_2.swf

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